



FOOD TECHNOLOGY FACT SHEET

Adding Value to OKLAHOMA

405-744-6071 • www.fapc.biz • fapc@okstate.edu

Edible Oil Quality

Nurhan Turgut Dunford

FAPC Oil/oilseed Specialist

Introduction

Majority of the edible oils used for cooking, frying and food formulations are derived from plant sources, specifically from oilseeds such as soybean, canola, sunflower seeds, cottonseed and peanuts. Edible vegetable oils are liquid at room temperature and comprised of mainly triacylglycerides that are made up of three fatty acids attached to a glycerol molecule through ester bonds (see fact sheet FAPC-196 Lipid Glossary). Physical, chemical and nutritional properties of vegetable oils vary significantly depending on the type of fatty acids present in the oil. Unsaturated fatty acid content of vegetable oils is significantly higher than that of animal fats. Although high unsaturated fatty acid content makes the oil healthier, highly unsaturated oils are prone to fast oxidation and quality deterioration during processing, handling and storage. Vegetable oils also contain other compounds in minor amounts that affect their quality and nutritional value. Phytosterols, tocopherols, waxes are some of the minor components of vegetable oils (less than 1 percent of the oil). These compounds are not regularly analyzed for oil quality evaluation purposes partly due to the complexity of the analytical protocols and the need for expensive instruments for testing.

Oil quality deteriorates by hydrolysis, oxidation and polymerization of the oil. Hydrolysis increases the amount of free fatty acids (FFA), mono- and diacylglycerols and glycerols in oils. Oxidation produces hydroperoxides and low molecular weight volatile compounds such as aldehydes, ketones, carboxylic acids, and short chain alkanes and alkenes. Dimers and polymers also are formed when oil is exposed to high temperatures during cooking and frying.

Monitoring and maintaining edible oil quality are paramount importance to ensure safety of the product for consumption. Although there is no official standard set for evaluating edible oil quality, FFA content, peroxide value (PV) and p-anisidine value (AV) are commonly used in industry to report edible oil quality. This fact sheet summarizes edible oil quality parameters used in industry and for research purposes.

Quality Parameters

Free Fatty Acid Content

Free fatty acids (FFA) are formed by the hydrolysis of oils (triacylglycerides). They are not bound or esterified to a glycerol molecule. Crude oils and fats in natural form, not refined, contain small amounts of FFA, which are usually removed during the refining process. FFA are not desirable in edible oils because when oils with high FFA content are used in foods, they lower the oxidative stability of the product, increase acidity and lead to off-flavor formation. The standard methods for determining FFA content of oils are American Oil Chemists' Society (AOCS) Ca-5a-40 [1] and European Union Regulation No.2568/91 [2]. There are other analytical methods, which are faster and suitable for small sample sizes [3,4]. Chemical analysis kits and various automated hand-held instruments also are available for analyzing FFA concentrations in oils [5].

Voluntary industry standard for FFA content in refined edible oil is ≤ 0.05 percent (based on oil weight). In the food industry, frying oils with FFA content exceeding 2 percent are either discarded or fresh oil is added to bring the FFA content down.

7. Dobarganes MC, Velasco J, Dieffenbacher A (2000) Determination of polar compounds, polymerized and oxidized triacylglycerols, and diacylglycerols in oils and fats: results of collaborative studies and the standardized method (Technical report). Pure and Applied Chemistry, vol 72. doi:10.1351/pac200072081563
8. Tarladgis B, Pearson AM, Dugan LR, Jr. (1962) The chemistry of the 2-thiobarbituric acid test for the determination of oxidative rancidity in foods. I. Some important side reactions. J Am Oil Chem Soc 39 (1):34-39. doi:10.1007/BF02633347
9. Tarladgis BG, Pearson AM, Jun LRD (1964) Chemistry of the 2-thiobarbituric acid test for determination of oxidative rancidity in foods. II.—formation of the tba-malonaldehyde complex without acid-heat treatment. Journal of the Science of Food and Agriculture 15 (9):602-607. doi:10.1002/jsfa.2740150904
10. St. Angelo A, Ory R, Brown L (1975) Comparison of methods for determining peroxidation in processed whole peanut products. Journal of the American Oil Chemists' Society 52 (2):34-35. doi:10.1007/BF02901817

The Oklahoma Cooperative Extension Service Bringing the University to You!

The Cooperative Extension Service is the largest, most successful informal educational organization in the world. It is a nationwide system funded and guided by a partnership of federal, state, and local governments that delivers information to help people help themselves through the land-grant university system.

Extension carries out programs in the broad categories of agriculture, natural resources and environment; home economics; 4-H and other youth; and community resource development. Extension staff members live and work among the people they serve to help stimulate and educate Americans to plan ahead and cope with their problems.

Some characteristics of Cooperative Extension are:

- The federal, state, and local governments cooperatively share in its financial support and program direction.
- It is administered by the land-grant university as designated by the state legislature through an Extension director.
- Extension programs are nonpolitical, objective, and based on factual information.

- It provides practical, problem-oriented education for people of all ages. It is designated to take the knowledge of the university to those persons who do not or cannot participate in the formal classroom instruction of the university.
- It utilizes research from university, government, and other sources to help people make their own decisions.
- More than a million volunteers help multiply the impact of the Extension professional staff.
- It dispenses no funds to the public.
- It is not a regulatory agency, but it does inform people of regulations and of their options in meeting them.
- Local programs are developed and carried out in full recognition of national problems and goals.
- The Extension staff educates people through personal contacts, meetings, demonstrations, and the mass media.
- Extension has the built-in flexibility to adjust its programs and subject matter to meet new needs. Activities shift from year to year as citizen groups and Extension workers close to the problems advise changes.

Oklahoma State University, in compliance with Title VI and VII of the Civil Rights Act of 1964, Executive Order 11246 as amended, and Title IX of the Education Amendments of 1972 (Higher Education Act), the Americans with Disabilities Act of 1990, and other federal and state laws and regulations, does not discriminate on the basis of race, color, national origin, genetic information, sex, age, sexual orientation, gender identity, religion, disability, or status as a veteran, in any of its policies, practices or procedures. This provision includes, but is not limited to admissions, employment, financial aid, and educational services. The Director of Equal Opportunity, 408 Whitehurst, OSU, Stillwater, OK 74078-1035; Phone 405-744-5371; email: eeo@okstate.edu has been designated to handle inquiries regarding non-discrimination policies; Director of Equal Opportunity. Any person (student, faculty, or staff) who believes that discriminatory practices have been engaged in based on gender may discuss his or her concerns and file informal or formal complaints of possible violations of Title IX with OSU's Title IX Coordinator 405-744-9154.

Issued in furtherance of Cooperative Extension work, acts of May 8 and June 30, 1914, in cooperation with the U.S. Department of Agriculture, Director of Oklahoma Cooperative Extension Service, Oklahoma State University, Stillwater, Oklahoma. This publication is printed and issued by Oklahoma State University as authorized by the Vice President, Dean, and Director of the Division of Agricultural Sciences and Natural Resources and has been prepared and distributed at a cost of 74 cents per copy. 0116 MG

Acid value

Acid Value is an important indicator of vegetable oil quality. Acid value is expressed as the amount of potassium hydroxide (KOH, in milligrams) necessary to neutralize free fatty acids contained in 1 g of oil. The standard method for determining acid value is AOCS Cd 3d-63 [1].

Peroxide Value

Peroxide value (PV) is an index used to quantify the amount of hydroperoxides present in fats and oils. Hydroperoxides, which are shown to be toxic to humans, are the primary oil oxidation products formed during the initial stages of oxidation. The standard method for PV determination is AOCS Cd 8-53 [1]. Although Fourier transform infrared (FTIR) and near-infrared (FT-NIR) spectroscopy methods developed for PV measurement have the advantages of analytical speed and automation, the instruments used for these tests are quite expensive and require extensive calibration [6]. Evaluating oil quality based on only PV can be misleading. Because low PV does not necessarily indicate low level of oxidation, it could be due to the advanced level of oil oxidation during which primary oxidation products are converted to secondary oxidation product (see the section for anisidine value) lowering PV but increasing AV of the oil. Hence, both PV and AV should be used for oil quality evaluation (see totox value).

Refined oils usually have PV of <1 meq/kg. Oils are considered oxidized when PV > 3 meq/kg.

p-Anisidine Value

p-Anisidine value (AV) is a measure of the secondary oxidation products that are formed by breakdown of the primary oxidation products during extensive oxidation. The secondary oxidation products are mainly aldehydes such as 2,4-dienals and 2-alkenals. AV is strongly correlated with overall oil odor intensity. The standard method for AV content determination in oils is the AOCS Cd 18-90 [1].

Refined oils should have AV of <5.

Polar Compounds

Presence of polar compounds in oil is one of the best indicators of heated oil quality. Polar compounds consist of dimeric and higher polymeric triacylglycerides formed through thermal polymerization of triacylglycerides, monomeric oxidized products, mono- and diacylglycerides and FFA formed through hydrolytic cleavage of triacylglycerides. The analysis of the polar

compounds is conducted by high-performance size exclusion chromatography, which allows the separation and quantification of polymeric compounds, dimers, oxidized triacylglycerides, mono- and diacylglycerides and FFA [7]. The standard method for analyzing total polar content (TPC) in frying oils is the AOCS Cd 20-91 [1]. As the name implies, this method determines the total polar compounds in the oil, not the individual polar compounds.

Regulations in some countries specify that oils should contain less than 25–27 percent TPC.

TBARS

This is one of the oldest tests used for evaluating lipid oxidation in foods. When heated under acidic conditions, thiobarbituric acid (TBA) reacts with a number of compounds including nucleic acids, amino acids, proteins, phospholipids and aldehydes to produce a pink chromophore that can be measured by UV or fluorescence detection. These substances are termed TBARS (thiobarbituric acid reacting substances). The extent of lipid oxidation is reported as TBA value, which corresponds to milligram of malonaldehyde equivalents per kilogram of sample or micromoles of malonaldehyde per gram of sample [8,9]. It is recommended other tests such as PV and AV also should be used when evaluating lipid quality in complex food systems due to the limitations of TBARS method. The standard method for analyzing TBA value is the AOCS Cd 19-90 [1].

Conjugated Dienes

Oxidation of polyunsaturated fatty acids results in an increase in the ultraviolet absorption of the product due to conjugate formation. Measurement of the content of conjugated dienes at 234 nm and conjugated trienes at 268 nm is a quick physical method, which may be helpful to assess the oxidative stability of vegetable oils [10].

Unsaponifiable Matter

Unsaponifiable matter (USM) fraction of vegetable oils naturally contains hydrocarbons, terpene alcohols, sterols, tocopherols and other phenolic compounds, which may act as oxidation inhibitors. Vegetable oils typically contain 0.5-2.5 percent USM while some others have higher amounts, 5-6 percent. USM of edible oils is used for their characterization and authentication of the products. The standard method for analyzing USM is the AOCS Ca 6b-53 [1].

Phospholipids

Phospholipid is a common name for lipids containing phosphoric acid or other phosphorus-containing acids in ester form such as glycerophospholipids (e.g. phosphatidic acid, phosphatidylcholine, phosphatidylethanolamine) or sphingophospholipids (e.g. sphingomyelin). Although these compounds (also called gum because of their gummy consistency in oil) have some health benefits and surfactant/emulsifier properties, they need to be separated from crude oil during the refining process, which is referred to as degumming. Otherwise, they impart a cloudy appearance and precipitate out of the oil during storage creating a unpleasant solid residue at the bottom of the containers and adversely affect the functionality of refined oils, i.e. cause foaming during frying. The phospholipid content of oils is commonly measured as phosphorous (AOCS Ca 19-86), which can be converted to phospholipids by using conversion factors calculated by using the phospholipid composition and the molecular weight of individual phospholipids present in the oil. The standard method for analyzing phospholipids is the AOCS Ca 12b-92 [1].

Refined oils have about 30 mg/kg phosphorous, while super degummed oils contain less than 10 mg/kg phosphorous.

Color

The color of crude and refined vegetable oils is an important factor in the determination of their market value. Removal of color pigments, which are extracted along with the oil from the seeds during the extraction process, is achieved during the oil refining referred to as bleaching. The color compounds in the oil mainly consist of carotenoids, chlorophyll, gossypol and related compounds. Chlorophyll is a sensitizer of photo-oxidation and promotes oil oxidation in the presence of light and decreases the oxidative stability of oils. Chlorophyll also acts as a catalyst poison during oil hydrogenation process. The color of the edible oils also can be an issue for food formulations adversely affecting the color of the final product which incorporated. Lovibond is the most common method (AOCS Wesson Cc 13b-45 and ISO Cc 13e-92) used to determine color of the commercial oils [1]. In Lovibond method color is expressed as red and yellow components. In general, fully refined oil may be 0.8 R (red) and 8.0 Y (yellow). Frying oils often are discarded when their Lovibond red color increases from 1.5-3.5 to 20-30.

Totox value

Totox value is used as an empirical assessment of oxidative deterioration based on PV and AV of an oil.

$$\text{Totox value} = 2 \times \text{PV} + \text{AV}$$

Conclusions

Commonly used edible oil industry voluntary quality parameters are FFA, PV and AV. National Oil Processors Association (NOPA) and American Oil Chemists' Society (AOCS) provide laboratory services that include AOCS methods, certified reference materials, quality reference materials and consulting. The AOCS Laboratory Proficiency Program (LPP), formerly the Smalley Check Sample Program, is the world's most extensive and respected collaborative proficiency testing program for oilseeds, oilseed meals, and edible oils and fats. More than 500 chemists who participate in this program use AOCS or similar standard methods for sample analyses and then compare their results with the test results generated at other laboratories using the same methods and samples and verify their quality control practices (<http://www.aocs.org/LabServices/content.cfm?ItemNumber=841>).

References

1. AOCS (2004) Official Methods and Recommended Practices. American Oil Chemists' Society, Champaign, IL.
2. Commission E Determination of the Free Fatty Acids in European Commission Regulation (EED) No. 2568/91. Official Journal of the European Communities L 248 (5.9.91):6
3. Mahesar SA, Sherazi STH, Khaskheli AR, Kandhro AA, Uddin S (2014) Analytical approaches for the assessment of free fatty acids in oils and fats. Anal Methods 6:4956-4963
4. Lowry RR, Tinsley IJ (1976) Rapid Colorimetric Determination of Free Fatty Acids. J Amer Oil Chem Soc 53 (4):470-472
5. Perkins EG, Erickson MD (1996) Deep Frying, Chemistry, Nutrition, and Practical Applications. AOCS Press, Champaign, IL
6. Li H, van de Voort FR, Ismail AA, Cox R (2000) Determination of peroxide value by fourier transform near-infrared spectroscopy. J Amer Oil Chem Soc 77 (2):137-142. doi:10.1007/s11746-000-0023-7